

Coated Filters Cat. No. 225-9035 Operating Instructions

Method & Chemical: IFV for toluene-2,4-diisocyanate and toluene-2,6-diisocyanate

Filter & Coating: 25-mm glass fiber filter coated with 1-(2-pyridyl)piperazine. Filters only, not loaded in cassettes

Prior to Sampling: Store at \leq 39.2 F (4 C). Limited shelf-life; check expiration date on packaging. **Ensure that filters are at room temperature before sampling.**

Sample Stability: Store samples at \leq 39.2 F (4 C) for up to 15 days.

For information on other available coated filters, go to www.skcinc.com.

SKC Limited Warranty and Return Policy

SKC products are subject to the SKC Limited Warranty and Return Policy, which provides SKC's sole liability and the buyer's exclusive remedy. To view the complete SKC Limited Warranty and Return Policy, go to skcinc.com/warranty.